
THUASNE

Lombacross Activity

fr	Ceinture de soutien lombaire renforcée	7
en	Reinforced lumbar support belt	8
de	Verstärkte LWS-Stabilisierungsothese	8
nl	Verstevigde lumbale rugbrace	9
it	Cintura di sostegno lombare rinforzata	10
es	Faja de sujeción lumbar reforzada	11
pt	Cinta de suporte lombar reforçada	12
da	Forstærket lændestøttebælte	13
fi	Vahvistettu lannerangan tukivyö	14
sv	Förstärkt stödbälte för ländryggen	15
el	Ενισχυμένη ζώνη οσφυϊκής υποστήριξης	16
cs	Vyztužený bederní pás	17
pl	Orteza lędźwiowa wzmocniona	18
lv	Stiprināta jostasvietas atbalsta josta	19
lt	Sutvirtintas palaikomasis juosmens diržas	20
et	Tugevdatud nimmepiirkonna tugivöö	20
sl	Ojačan oporni križni pas	21
sk	Vystužený podporný driekový pás	22
hu	Megerősített ágyéki fűző	23
bg	Подсилен лумбален колан	24
ro	Centură de susținere lombară întărită	25
ru	Усиленный опорный поясничный корсет	26
hr	Ojačani lumbalni potporni pojas	27
zh	加强型腰部支撑带	28
ar	حزام مقوى لدعم أسفل الظهر	29

	52 - 62 cm	1
	63 - 74 cm	2
	75 - 89 cm	3
	90 - 105 cm	4
	106 - 120 cm	5
	121 - 140 cm	6

fr	Tour de taille	Hauteurs : 21 cm - 26 cm
en	Waist measurement	Heights: 21 cm - 26 cm
de	Taillenumfang	Höhen: 21 cm - 26 cm
nl	Tailleomvang	Hoogtes: 21 cm - 26 cm
it	Circonferenza vita	Altezze: 21 cm - 26 cm
es	Contorno de cintura	Alturas: 21 cm - 26 cm
pt	Perímetro de cintura	Alturas: 21 cm - 26 cm
da	Taljemål	Højde: 21 cm - 26 cm
fi	Vyötärönympärys	Korkeudet: 21-26 cm
sv	Midjemått	Höjd: 21 cm - 26 cm
el	Περίμετρος μέσης	Ύψη: 21 εκ. - 26 εκ.
cs	Obvod pasu	Výšky: 21 cm - 26 cm
pl	Obwód w talii	Wysokości: 21 cm-26 cm
lv	Vidukļa apkārtmērs	Augstumi: 21 cm - 26 cm
lt	Juosmens apimtis	Aukštis: 21-26 cm
et	Vööümbmõõt	Kõrgused: 21 cm - 26 cm
sl	Obseg pasu	Višina: 21 cm-26 cm
sk	Obvod pása	Výšky: 21 cm 26 cm
hu	Derékbőség	Magasság: 21 cm - 26 cm
bg	Обиколка на талията	Височини: 21 cm и 26 cm
ro	Circumferința taliei	Înălțimi: 21 cm - 26 cm
ru	Обхват талии	Высота: 21-26 см
hr	Opseg struka	Visine: 21 cm do 26 cm
zh	腰圍	高度：21 cm - 26 cm
ar	محيط الخصر	الارتفاعات: 21 سم - 26 سم

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REINFORCED LUMBAR SUPPORT BELT

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Aerated and breathable Combitec fabric. ❶

Optimized back comfort due to a fleecy-foam complex.

Abdominal panels with self-fastening surface and finger-loop. ❷

Additional strap with thumb-loops. ❸

Conformable dorsal stays. ❹

Textile components: Polyamide - Polyester - Elastane
Non-textile components: Brushed tempered steel - Polyoxymethylene - Polypropylene - PVC

Properties/Mode of action

The device relieves low back pain, supporting the lumbar spine: compartmental lumbar unloading is therefore achieved via increased intra-abdominal pressure.

Indications

Non-specific low back pain (acute, subacute and chronic phases).

Spinal stenosis.

Contraindications

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Precautions

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

Do not wear the product in a medical imaging machine.

Do not wear the product while sleeping.

Consult a health professional if discomfort or pain occurs.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Undesirable side-effect

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Position the belt around the waist, placing your hands in the finger loop located on the front panel. ❶

The belt must be centered on the spine: the dorsal stays should be positioned on either side of the spine.

The bottom of the belt must stop at the top of the crease between the buttocks.

The inside label must be positioned upwards.

Hold the left end flat using your left hand.

Pull the finger loop on the abdominal panel simultaneously and evenly. ❷

Fasten the belt, right side over left side, adjusting to the desired compression. ❸

For an adapted compression, you can use the additional strap.

Pull the thumb loops on the additional strap simultaneously and evenly. ❹

Attach the self-fastening parts onto the centre of the abdominal panels of the belt. ❺

The level of compression can be adjusted as required throughout the day.

Care

Close the self-fastening tabs before washing. Machine washable at 40 °C. Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2020301 (2021-02)
Thuasne SAS – SIREN/RCS Nanterre 542 091 186 – capital 1 950 000 euros