
THUASNE

Lombacross Activity

| | | |
|----|---|----|
| fr | Ceinture de soutien lombaire renforcée avec compression homogène..... | 5 |
| en | Reinforced lumbar support with homogeneous compression | 6 |
| de | Verstärkte LWS-Stabilisierungsorthese mit gleichmäßiger Kompression..... | 7 |
| nl | Verstevigde lumbale rugbrace met gelijkmatige compressie..... | 8 |
| it | Cintura di sostegno lombare rinforzata a compressione omogenea | 9 |
| es | Faja de sujeción lumbar reforzada con compresión homogénea..... | 10 |
| pt | Cinta de suporte lombar reforçada com compressão homogénea..... | 11 |
| da | Forstærket lændestøttebælte med ensartet kompression..... | 12 |
| fi | Vahvistettu lannerangan tukivyö - tasainen puristus..... | 13 |
| sv | Förstärkt ryggörcdel för ländryggen med jämn kompression..... | 14 |
| el | Ζώνη ενισχυμένης οσφυϊκής υποστήριξης με ομοιόμορφη συμπίεση | 15 |
| cs | Vyztužený bederní pás s rovnoměrnou kompresí..... | 17 |
| pl | Wzmocniona orteza lędźwiowa z równomiernym naciskiem | 18 |
| lv | Stiprināta muguras jostas daļas atbalsta josta ar vienmērīgu kompresiju | 19 |
| lt | Sutvirtintas palaikomasis vienodo suspaudimo juosmens diržas | 20 |
| et | Nimpepiirkonna tugevdatud, ühtlase kompressiooniga tugivöö | 21 |
| sl | Ojačan oporni križni pas s homogeno kompresijo | 22 |
| sk | Spevnený driekový pás s rovnomernou kompresiou | 23 |
| hu | Megerősített ágyéki fűző homogén kompresszióval | 24 |
| bg | Подсилен лумбален колан с хомогенна компресия..... | 25 |
| ro | Cintură de susținere lombară întărită cu compresie omogenă..... | 26 |
| ru | Усиленный опорный поясничный корсет с равномерной компрессией..... | 27 |
| hr | Ojačani lumbalni potporni pojasi s homogenom kompresijom | 28 |
| zh | 均匀施压的加强型腰部支撑带 | 29 |
| ar | حزام مقوّي لدعم أسفل الظهر بميزة الضغط المتجانس | 31 |

| | | |
|---|---|--------------|
|  | 1 | 60 - 80 cm |
| | 2 | 75 - 95 cm |
| | 3 | 90 - 110 cm |
| | 4 | 105 - 125 cm |
| | 5 | 120 - 140 cm |








| | | |
|-----------|----------------------|---------------------------|
| fr | Tour de taille | Hauteurs : 21 cm - 26 cm |
| en | Waist circumference | Heights: 21 cm - 26 cm |
| de | Taillenumfang | Höhen: 21 cm - 26 cm |
| nl | Tailleomvang | Hoogtes: 21 cm - 26 cm |
| it | Circonferenza vita | Altezze: 21 cm - 26 cm |
| es | Contorno de cintura | Alturas: 21 cm - 26 cm |
| pt | Perímetro de cintura | Alturas: 21cm - 26cm |
| da | Taljemål | Højde: 21cm - 26cm |
| fi | Vyötärönympärys | Korkeudet: 21 cm - 26 cm |
| sv | Midjemått | Höjd: 21cm - 26cm |
| el | Περίμετρος μέσης | Ύψη: 21 εκ. - 26 εκ. |
| cs | Obvod pasu | Výšky: 21 cm - 26 cm |
| pl | Obwód w talii | Wysokości: 21 cm - 26 cm |
| lv | Vidukļa apkārtmērs | Augstumi: 21 cm - 26 cm |
| lt | Juosmens apimtis | Aukštis: 21 cm - 26 cm |
| et | Vööümbermõõt | Kõrgused: 21 cm - 26 cm |
| sl | Obseg pasu | Višina: 21 cm - 26 cm |
| sk | Obvod pása | Výšky: 21 cm - 26 cm |
| hu | Derékбőсég | Magasság: 21 cm - 26 cm |
| bg | Обиколка на талията | Височини: 21 cm - 26 cm |
| ro | Circumferința taliei | Înălțimi: 21 cm - 26 cm |
| ru | Обхват талии | Высота: 21 см - 26 см |
| hr | Opseg struka | Visine: 21 cm - 26 cm |
| zh | 腰围 | 高度 : 21 cm -26 cm |
| ar | محيط الخصر | الارتفاعات: 21 سم - 26 سم |

| | | | |
|-----------|--------------|----------------|--------------------------------|
| fr | Compression | Maintien | Education posturale |
| en | Compression | Support | Postural education |
| de | Kompression | Halt | Erlernen der Körperhaltung |
| nl | Compressie | Fixatie | Voorlichting over houding |
| it | Compressione | Sostegno | Educazione posturale |
| es | Compresión | Sujeción | Educación postural |
| pt | Compressão | Manutenção | Educação postural |
| da | Kompression | Støtte | Postural genoptræning |
| fi | Kompressio | Tuki | Asentokuntoutus |
| sv | Kompression | Stöd | Postural träning |
| el | Συμπίεση | Συγκράτηση | Σωστή στάση σώματος |
| cs | Komprese | Stabilizace | Správné držení těla |
| pl | Kompresja | Podtrzymywanie | Edukacja dotycząca postawy |
| lv | Kompresija | Atbalsts | Stājas korekcija |
| lt | Kompresija | Palaikymas | Laikysenos taisymas |
| et | Kompressioon | Hooldus | Rühitreening |
| sl | Kompresija | Oprjem | Učenje telesne drže |
| sk | Kompresia | Podpora | Cvičenie posturálnych svalov |
| hu | Szorítás | Támasz | Tartásjavító oktatás |
| bg | Компресия | Поддръжка | Приучаване към правилна стойка |
| ro | Compresie | Susținere | Educație posturală |
| ru | Компрессия | Стабилизация | Исправление осанки |
| hr | Kompresija | Potpora | Učenje o držanju |
| zh | 压紧 | 稳定 | 姿势联系 |
| ar | الضغط | الدعم | تدريب وضعي |

fr

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en **REINFORCED LUMBAR SUPPORT** **WITH HOMOGENEOUS** **COMPRESSION**

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Aerated and breathable Combitex fabric. ①

Abdominal panels with self-fastening surface and finger-loops. ②

Additional strap with finger-loops. ③

Comfortable dorsal stays. ④

Composition

Textile components: polyamide - polyester - elastane - polyurethane.

Non-textile components: brushed tempered steel - polyoxymethylene - polypropylene - polyethylene.

Properties/Mode of action

The device relieves low back pain by supporting the spine through an increased intra-abdominal pressure which eases the load on the intervertebral discs.

Indications

Non-specific low back pain (acute, subacute and chronic phases). Prevention of recurrent low back pain in combination with return to activity and/or therapeutic information.

Spinal stenosis.

Spine stability.

Resumption of activity.

Contraindications

Do not use in the event of known allergy to any of the components.

*Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Precautions

Verify the product's integrity before every use.

Do not use the product if it is damaged.

Wear the device above thin clothing.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not wear the product while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Measure the waist circumference of the patient at the mid-distance between the iliac crests and floating ribs.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not wear the product in a medical imaging machine.

Undesirable side-effect

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time. ⑥

Position the belt around the waist, using the finger loops located on the front panels. ⑦

The Thuasne logo indicates which way up the belt should be put on.

Pull the finger loops on the abdominal panels simultaneously and evenly.

Fasten the belt, right side over left side, adjusting to the desired compression. ⑧

The belt must be centred on the spine: the dorsal stays should be positioned on either side of the spine. ⑨

The bottom of the belt must stop at the top of the crease between the buttocks. ⑩

For a reinforced compression, tighten the additional strap.

Pull the finger loops on the additional strap simultaneously and evenly. ⑪

Attach the self-fastening parts of the additional strap onto the belt's abdominal panels or place them on top of each other. ⑫

The level of compression can be adjusted as required throughout the day.

Care/Maintenance

Product can be washed in accordance with the instructions shown on this leaflet and on the label.

Close the self-fastening tabs before washing.

Machine washable at 30°C (delicate programme). The device has been tested over 15 wash cycles (average: one wash per month). Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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Suosittelimme, että terveydenhuollon ammattilainen säätää pystytuet potilaan selälle sopiviksi ennen ensimmäistä käyttökertaa. ©
THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091 186